

# **DEXUP Mechanical CPR Device (MCPR) – FAQ**

## **1. General Overview**

### **Q1. What is the DEXUP Mechanical CPR (MCPR) device?**

The DEXUP MCPR is an automated CPR solution that delivers reliable, high-quality chest compressions during cardiac arrest, improving effectiveness and reducing manual effort.

### **Q2. What is the primary clinical objective of this device?**

The device is designed to:

- Deliver uniform and guideline-compliant chest compressions
- Minimize interruptions during CPR
- Reduce physical strain on healthcare providers
- Improve overall resuscitation outcomes

### **Q3. In what environments can the DEXUP MCPR device be deployed?**

The device is suitable for a wide range of emergency response environments, including:

- Hospitals (ICU, Emergency Departments, Operation Theaters)
- Ambulances and pre-hospital emergency care units
- Corporate offices and large workplaces with occupational health programs
- Public infrastructure such as airports, railway stations, and metro stations
- Industrial facilities and high-risk work zones
- Disaster response and mass casualty management scenarios

Its portability and rapid deployment capability make it highly effective in both medical and non-medical emergency settings.

## **2. Features & Functional Capabilities**

### **Q4. What are the key operational features of the device?**

- Automated, high-quality chest compressions
- One-button rapid activation system
- Compact and ergonomically optimized design
- Intelligent chest height detection
- Multiple compression modes
- Controlled and gradual compression initiation

- Simplified and intuitive user interface

#### **Q5. How does the device ensure consistency in CPR delivery?**

The system uses automated mechanical actuation to maintain:

- Consistent compression depth
- Stable compression rate
- Reduced variability compared to manual CPR

#### **Q6. What is meant by “automated chest height detection”?**

This feature enables the device to automatically assess the patient’s chest dimensions and adjust compression depth accordingly, ensuring optimal performance across different body types within the supported range.

#### **Q7. What compression modes are supported?**

The device supports:

- Continuous compression mode
- 30:2 compression-to-ventilation ratio
- 15:2 compression-to-ventilation ratio

These modes align with standard resuscitation protocols.

#### **Q8. What is the benefit of gradual compression initiation?**

Gradual initiation ensures that compressions begin progressively rather than abruptly, helping:

- Reduce the risk of injury
- Improve mechanical adaptation to the patient’s chest

### **3. Technical Specifications & Performance**

#### **Q9. What is the model and classification of the device?**

- Model: **DEXUP MCPR**
- Device Type: Mechanical CPR System

#### **Q10. What are the electrical input requirements?**

- Operating Voltage: 100–240V AC
- Frequency: 50/60 Hz
- Power Consumption: 180 VA

### **Q11. What are the battery performance parameters?**

- Battery Voltage: 24V
- Capacity: 4500 mAh
- Operating Time: Up to 120 minutes (continuous use)

### **Q12. What are the charger output specifications?**

- Output Voltage: 24V
- Output Current: 3A

### **Q13. What are the device dimensions and weight?**

- Dimensions: 65 cm (Height) × 54 cm (Width)
- Net Weight: 8.5 kg

This ensures portability without compromising structural stability.

### **Q14. What level of environmental protection does the device offer?**

- IP Rating: IPX3 (excluding charger)

This indicates protection against light water spray, making it suitable for emergency environments.

## **4. Patient Applicability**

### **Q15. What patient size range is supported by the device?**

- Sternum Height: 17.5 cm to 29 cm
- Chest Width: Up to 48 cm

### **Q16. Is the device suitable for pediatric use?**

No. The device is not recommended for pediatric patients due to dimensional and safety limitations.

### **Q17. Can the device be used for bariatric (obese) patients?**

Use is not recommended in cases of extreme obesity due to potential challenges in positioning and compression effectiveness.

## **5. Contraindications & Safety Considerations**

**Q18. What are the major contraindications for using the device?**

The device should not be used in patients with:

- Severe chest trauma or deformities
- Flail chest
- Open chest wounds
- Significant thoracic injuries

**Q19. Is the device safe for post-surgical cardiac or thoracic patients?**

No. Use may lead to:

- Surgical wound disruption
- Bleeding complications
- Compromise of surgical repairs

**Q20. Can it be used in patients with implanted cardiac or thoracic devices?**

Caution is advised due to the risk of:

- Device displacement
- Functional interference

Clinical evaluation is recommended prior to use.

**Q21. What precautions are required for pregnant patients?**

Use requires modified positioning due to the risk of inferior vena cava compression. Clinical judgment is essential.

## **. Clinical & Operational Advantages**

**Q22. What are the primary advantages of the DEXUP MCPR device?**

- Ensures consistent, guideline-compliant compressions
- Significantly reduces rescuer fatigue
- Enhances operational efficiency during resuscitation
- Enables continuous CPR during patient transport
- Improves overall workflow in emergency scenarios

**Q23. How does the device improve clinical outcomes?**

By maintaining uninterrupted and high-quality compressions, the device supports:

- Improved blood circulation during cardiac arrest
- Reduced human error
- Better adherence to CPR protocols

## **7. Usability & Deployment**

### **Q24. How user-friendly is the device?**

The device is designed for rapid adoption with:

- One-button operation
- Automated calibration
- Minimal training requirements

### **Q25. How quickly can the device be deployed in an emergency?**

The system supports rapid deployment, enabling immediate initiation of CPR within seconds.

## **. Portability & Durability**

### **Q26. Is the device suitable for mobile and field use?**

Yes. Its lightweight and compact construction make it ideal for:

- Ambulances
- Field emergencies
- Public access emergency systems

### **Q27. Is the device resistant to environmental exposure?**

With an IPX3 rating, it is protected against light splashes, though it is not designed for heavy water exposure or submersion.

## **9. Manufacturer & Support**

### **Q28. Who manufactures the device?**

The device is developed and supplied by **Codex Healthcare**.

### **Q29. How can support or sales inquiries be made?**

- Website: [www.dexup.in](http://www.dexup.in)

- Phone: +91 73853 04200
- Email: [sales@dexup.in](mailto:sales@dexup.in)